Digital Bathroom Scale Model HD-HAILI-007

EN INSTRUCTIONS FOR USE

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


EN INSTRUCTIONS FOR USE

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Representative Tanita Europe B.V.
Hoogoordreef 56-E
1101 BE Amsterdam
The Netherlands
Tel: +31(0)20 560 2970
Fax: +31(0)20 560 2988
www.tanita.eu

Tanita Corporation of America Inc.
2625 South Clearbrook Drive
Arlington Heights, IL 60005 USA
Tel: +1-847-640-9241
Fax: +1-847-640-9261
www.tanita.com

Tanita Health Equipment H.K. Ltd.
Unit 301-303 3/F Wing On Plaza,
62 Mody Road, Tsimshatsui East,
Kowloon, Hong Kong
Tel: +852-2834-3917
Fax: +852-2838-8667
www.tanita.asia

ANITA India Private Limited
Level 8, Vibgyor Towers, G Block,
C-62
Bandra Kurla Complex, Bandra
(East),
Mumbai 400051, INDIA
Tel: +91-22-4090-7174
Fax: +91-22-4090-7178

Manufacturer Tanita Corporation
1-14-2, Maeno-Cho, Itabashi-ku,
Tokyo, Japan 174-8630
Tel: +81(0)-3-3968-2123
Fax: +81(0)-3-3967-3766
www.tanita.co.jp
ISO 9001 Certified

EN INSTRUCTIONS FOR USE

1. Specifications

Capacity: 150kg
Graduation: 0,1kg
Battery Type: AA (R6) x 2

2. Parts (Fig. 1)

a) Body b) Display c) Battery (Bottom)

3. Activating The Scale (Only When Installing New Batteries)

- 1) Insert batteries. (Remove the insulation sheet inserted at shipment from the battery contacts.)
- 2) Place the scale on a firm even surface, making sure all four feet make contact at the same time.
- 3) With your foot tap the middle of the platform to activate the scale. On the LCD display "88888" will appear - wait until it changes to "0.0" before stepping onto the platform.
 - * If "Err" appears on the display, redo the steps from 2).

4. Weighing

- 1) Step onto the platform and stand still.
 - * When weighing items of 12kg or less, first load a weight of 12kg or more and then weigh the item after "0.0kg" is displayed.
- 2) When the weight is determined, the measured value flashes on the display.
- 3) The scale will turn off automatically after 20 seconds.
 - * The scale should be used on a firm surface for accurate weighing, not recommended for use on carpeted surfaces.
 - * The LCD will display "OL", if the weight is over 150 kg and will turn off after 30 seconds.
 - * Whenever moving the scale to another location make sure that it is turned off.
 - * Do not move the unit or place anything on the unit when "OFF" is displayed.

5. Replacing the Batteries

When the low battery sign (Lo) appears on the LCD, remove the current batteries and replace them with new set as indicated in Fig 2.

6. Care Instructions

- 1) Be careful not to open the scale or tamper with its mechanism or circuitry. These actions will void all warranty coverage.
- 2) Take care not to spill water or chemical liquids onto the scale. Wipe the scale with a damp cloth periodically.
- 3) When not used for long periods of time, remove the batteries to prevent damage.
- 4) Do not use this scale while operating transmitters (such as a cell phone), as this may affect your result.
- 5) Do not step on the platform when wet.
- 6) Do not drop any objects onto the platform.

- 7) If the scale is placed against or under something, it may not show the correct weight when it is used again. Should this happen,
 1. Place the scale on a firm even surface.
 2. With your foot tap the middle of the platform to activate the scale.
 3. Once again, with your foot tap the middle of the platform until the display will disappear.
 4. Step on to the platform.

Federal Communications Commission and Canadian ICES Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and Canadian ICES-003. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular in tallation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures.

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio or television technician for help.

Modifications

The FCC requires the user to be notified that any changes or modifications made to this device that are not expressly approved by Tanita Corporation may void the user's authority to operate the equipment.

This device meets the EMC Directive 2004/108/EC and the RoHS Directive 2011/65/EU.

Notice:
These marks can be found out on contaminant-containing Batteries:
Pb Pb = Battery contains lead,
Cd Cd = Battery contains cadmium,
Hg Hg = Battery contains mercury.

⚠ Not allowed to mix batteries with consumer wastes!
As consumer you are legally bound to return used or discharged batteries. You can deposit your old batteries at the public collecting points in your town, or wherever the corresponding batteries are sold and specifically marked collecting boxes have been set up. In case of scrapping the apparatus, the batteries should be removed from it and deposited at the collecting points as well.

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